

ULTRABLOX® X-Ray Attenuating Cream



ULTRABLOX® Attenuating Cream allows you to provide the necessary care to your patients while significantly reducing dose exposure to your own hands. The cream protects the hands while retaining dexterity and tactile feel, providing a practical and efficient alternative to heavy lead gloves.



- Patented formula reduces scatter radiation exposure by up to 85%.
- Provides increased clinician safety through effective radiation dose reduction.
- Retains dexterity of the fingers and preserves tactile feel.
- Smooth and easy to use; can be applied directly to hands or encapsulated between two sets of gloves.
- Can be easily incorporated into physician's prep procedure.
- Biocompatible and non-toxic.
- Each 50 gram tube of sterile attenuating cream is contained inside a peel pouch.
- Packaged 15 tubes per box.



ULTRABLOX Cream has received a 510(k) from the FDA. U.S. Patent No. 9,114,121.



Ordering Information for BLOXR® Products

•ULTRABLOX® Attenuating Cream

P/N	Qty.
100036	15/box

•Machine Washable, Non-Toxic XPF® Attenuating Shields

Available in: Black (-BLK), Blue (-BLU), Green (-GRN), Red (-RED), Purple (-PUR), Grey (-GRY)

CAP			THYROID COLLAR		
P/N	Size	Head	P/N	Size	Neck
100098	S	21" to 23"	100007	S	10" to 15"
100099	M	22" to 24"	100008	M	16" to 18"
100100	L	23" to 25"	100009	L	19" to 21"

SKIRT					VEST				
P/N	Size	Waist	Hips	Length	P/N	Size	Chest	Low Waist	Length
100106	XS	26" to 29"	36"	19"	100096	XS	30" to 32"	32" to 34"	21"
100027	S	29" to 32"	38"	20"	100030	S	32" to 35"	37" to 39"	22"
100028	M	32" to 34"	41"	21"	100031	M	36" to 41"	39" to 42"	24"
100029	L	34" to 36"	44"	22"	100032	L	41" to 46"	42" to 44"	25"
100039	XL	37" to 40"	49"	23"	100040	XL	46" to 51"	44" to 47"	27"
100102	XXL	40" to 43"	52"	24"	100097	XXL	48" to 53"	48" to 51"	28"

FRONTAL APRON					WRAP APRON				
P/N	Size	Chest	Length	Max Hgt.	P/N	Size	Chest	Length	Max Hgt.
100019	XS	30" to 34"	32"	5' 0"	100107	XS	30" to 32"	32"	5' 0"
100010	S	34" to 38"	35"	5' 5"	100103	S	32" to 35"	35"	5' 5"
100011	M	38" to 42"	40"	5' 10"	100104	M	36" to 41"	40"	5' 10"
100012	L	42" to 46"	44"	6' 3"	100105	L	41" to 46"	44"	6' 3"
100038	XL	46" to 50"	46"	6' 8"	100108	XL	46" to 51"	46"	6' 8"

DROP APRON					PATIENT PROTECTION SHIELDS			
P/N	Size	Chest	Length	Max Hgt.	P/N	Size	P/N	Size
100041	S	34" to 38"	37"	5' 5"	100201	18" x 18"	100210	24" x 24"
100042	M	38" to 42"	42"	5' 10"	100202	20 x 20"	100211	24" x 36"
100043	L	42" to 46"	47"	6' 3"	Contact us for a full list of available sizes			

•Disposable, Non-Toxic XPF® Attenuating Garments

Not washable; light blue

CAP				THYROID COLLAR			
P/N	Pkg.	Size	Head	P/N	Pkg.	Size	Neck
100004	12/box	S	21" to 23"	100001	12/box	S	10" to 15"
100005	12/box	M	22" to 24"	100002	12/box	M	16" to 18"
100006	12/box	L	23" to 25"	100003	12/box	L	19" to 21"

•Apron Accessories

HANGERS		ROLLING RACKS		PEG RACKS Wall Mount	
P/N	Type	P/N	Type	P/N	Type
7XHANG	Standard	7XSPID	Spider style	7XPEG02	2 Pegs
7XHANGBA	Brushed Aluminum	7XVRACK	Valet style	7XPEG04	4 Pegs
SWING RACK Wall Mount		7XVRACK10	Valet style with 10 hangers	7XPEG08	8 Pegs
7XSWNGR	Right				
7XSWNGL	Left				
7XSWNGF	Full Set				

•Radiation Protection Eyeglasses

Visit bloxr.com for an up-to-date list of available frames and options



BLOXR

BLOXR Solutions, LLC

P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 | Fax: 877-254-4888
order@bloxr.com | www.bloxr.com

U.S. Patent No. 8,754,389, 8,993,989 and 9,114,121.
©2015, 2016 BLOXR Solutions

79021A



BLOXR



- ULTRABLOX® Attenuating Cream
- Washable & Flexible XPF® Attenuating Shields

CE 0120

FDA 510(k) Clearance

Registered with Health Canada



**MACHINE WASHABLE XPF®
Radiation Attenuating Shields**



BLOXR

- Effective Protection
- Washable
- Lightweight
- Durable
- Bendable
- Non-Toxic
- No Metals



Optional embroidery is available.



CAP
Ties securely in back. Wicking material inside the cap and a mesh top provide maximum comfort. Provides 0.5mm lead equivalency around the circumference of the head.



THYROID COLLAR
Secures with hook and loop closures at the back of the neck for a comfortable fit. Provides 0.5mm lead equivalency over the thyroid area.



FRONTAL APRON
Equipped with a belt for a secure fit. Provides 0.5mm lead equivalent protection in the front.



*** NEW APRON ACCESSORIES**



VEST
Each garment provides 0.5mm lead equivalent protection in the front where the panels overlap and 0.25 Pb equivalent coverage on the back.



SKIRT

- BLOXR® XPF® garments provide 0.5mm Pb equivalency with no compromise in protection across the full kVp spectrum.
- Unique, patented bi-layer provides lead-equivalent protection against scatter radiation at almost half the weight of traditional lead, helping to alleviate back pain and fatigue.
- Can be bent, folded and flexed without cracking. You can hang them over a chair or sit down while wearing an XPF garment.
- Durable construction -- garments can be machine washed and wiped down for cleanliness.
- Since there are no lead or heavy metals used, there are no special requirements for disposal. Three year warranty.



WRAP APRON
Inner belt for secure fit. Provides 0.5mm lead equivalent protection in the front where the panels overlap and 0.25 Pb equivalent coverage on the back.

DROP APRON
Allows for quick removal. Provides 0.5mm lead equivalent protection in the front.



*** NEW PATIENT SHIELDS**
XPF® Patient Protection Shields provide 0.5mm lead equivalent protection. Place over a patient's radio-sensitive organs. Various lengths and widths available.



CUSTOM OPTIONS:
Do you require non-standard sizing? Contact your BLOXR representative for details.

Personalize with (A) embroidery, add (B) tethering, add a (C) support belt, include a (D) dosimeter pocket, protect your arm with a (E) sleeve, or add an anti-slip liner.



*** NEW - Radiation Protection Eyeglasses**

- Lenses provide 0.75 mm lead equivalency.
- Various frame styles and colors are available; optional Rx and side shields. Visit our website to view most current options.
- Eyeglasses carry a lifetime warranty on defects and workmanship.



Shown: Boston Blue, Crusier Tortoise, Nike Brazen Gray

DISPOSABLE XPF® Radiation Attenuating Caps & Collars

- Made of light blue polypropylene; not washable.
- Completely non-toxic; contains no lead or heavy metals.
- Packaged 12 per box.



Cap ties securely in back. Provides 0.5mm lead equivalency around the circumference of the head.



Thyroid collar secures with hook and loop closures at the back of the neck. Provides 0.5mm lead equivalency over the thyroid area.

ABOUT XPF® TECHNOLOGY



BLOXR

BLOXR® XPF® garments feature a patented, lightweight bi-layer constructed of barium sulfate and bismuth oxide. The two layers work together to provide the best attenuation by weight as compared to other common attenuating materials. XPF filters out the lower energy, more harmful part of the X-ray beam, thereby providing increased safety while reducing dose. XPF is non-toxic and is free of lead and heavy metals so garments can be disposed of without special consideration. BLOXR XPF has received 510(k) clearance from the FDA.

U.S. Patent No. 8,754,389 and 8,993,989.



XPF Garments are flexible and bendable.