



BLOXR

Featured Dental Products

For scatter radiation protection

XPF® Patient Poncho

Minimizes a patient's exposure to ionizing radiation while providing anatomical access for the clinician. Poncho is lightweight so it doesn't put pressure on the patient. Built in flap covers the thyroid area. Extenders can be added for extra length for taller patients. 0.5 mm Pb lead equivalency in front.

Patient Poncho		Extender Sets (optional)		
PN*	Size	Size	PN*	Lengths
100231	#1 Small Child	For #1	100240	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm)
100232	#2 Child	For #2	100241	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)
100233	#3 Small Adult	For #3	100242	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)
100234	#4 Adult	For #4	100243	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)
100235	#5 Large Adult	For #5	100244	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)



FRONT
shown with extender



BACK

XPF® Patient Shields

Designed for placement over radio-sensitive organs and gonads to protect the patient from harmful ionizing scatter radiation. 0.5 mm Pb lead equivalency in front.

PN*	Size	PN*	Size
100201	18" x 18" (46 x 46 cm)	100210	24" x 24" (61 x 61 cm)
100202	20" x 20" (51 x 51 cm)	100211	24" x 36" (61 x 91 cm)
100203	14" x 16" (36 x 41 cm)	100212	24" x 46" (61 x 117 cm)
100204	11" x 34" (28 x 86 cm)	100213	8" x 10" (20 x 25 cm)
100206	16" x 24" (41 x 61 cm)	100214	12" x 12" (31 x 31 cm)
100207	18" x 24" (46 x 61 cm)	100215	16" x 12" (41 x 31 cm)
100208	20" x 24" (51 x 61 cm)	100216	12" x 14" (31 x 36 cm)



XPF® Frontal Aprons

Provides 0.5 mm lead equivalent protection against scatter radiation in the front. Equipped with a belt, padded shoulders, one breast pocket and two hip pockets for convenience and comfort. Lightweight, machine washable.

PN*	Size	Chest	Length	Max Hgt
100019	XS	in: 30" - 34" cm: 76 - 86	in: 32" cm: 81	5' 0" cm: 152
100010	S	in: 34" to 38" cm: 86 - 97	in: 35" cm: 89	5' 5" cm: 165
100011	M	in: 38" - 42" cm: 97 - 107	in: 40" cm: 102	5' 10" cm: 178
100012	L	in: 42" - 46" cm: 107 - 117	in: 44" cm: 112	6' 3" cm: 191
100038	XL	in: 46" - 50" cm: 117 - 127	in: 46" cm: 117	6' 8" cm: 203
100045	XXL	in: 50" - 56" cm: 127 - 142	in: 47" cm: 119	6' 10" cm: 208



XPF® Thyroid Collars

Comfortably fits around the neck. Lead equivalency 0.5 mm Pb over thyroid. Machine washable. Secures with a hook and loop closure at the back of the neck. Patented XPF® Technology.



PN*	Size	Neck
100007	S	in: 10" - 14" cm: 25 - 35
100008	M	in: 15" - 17" cm: 38 - 43
100009	L	in: 18" - 20" cm: 46 - 51
100049	XL	in: 21" - 23" cm: 53 - 58



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com

*Add color code