



BLOXR

Featured Veterinary Products

For scatter radiation protection

XPF® Frontal Aprons

Provides 0.5 mm lead equivalent protection against scatter radiation in the front. Equipped with a belt, padded shoulders, one breast pocket and two hip pockets for convenience and comfort. Lightweight, helps reduce back pain and the limiting mobility of heavy lead aprons. Machine washable. Patented XPF® Technology.

PN	Size	Chest	Length	Max Hgt
100019	XS	in: 30" - 34" cm: 76 - 86	in: 32" cm: 81	5' 0" cm: 152
100010	S	in: 34" to 38" cm: 86 - 97	in: 35" cm: 89	5' 5" cm: 165
100011	M	in: 38" - 42" cm: 97 - 107	in: 40" cm: 102	5' 10" cm: 178
100012	L	in: 42" - 46" cm: 107 - 117	in: 44" cm: 112	6' 3" cm: 191
100038	XL	in: 46" - 50" cm: 117 - 127	in: 46" cm: 117	6' 8" cm: 203
100045	XXL	in: 50" - 56" cm: 127 - 142	in: 47" cm: 119	6' 10" cm: 208



XPF® Mitts

One size fits all.

CLOSED Mitts have a sewn in slit to allow fingers to be freed if desired.



Closed PN 100048 Style Closed, pair

OPEN Mitts enable easy movement of fingers.



Open PN 100037 Style Open, pair

XPF® Thyroid Collars

Comfortably fits around the neck. Lead equivalency 0.5 mm Pb over thyroid. Machine washable. Secures with a hook and loop closure at the back of the neck. Patented XPF® Technology.

PN	Size	Neck
100007	S	in: 10" - 14" cm: 25 - 35
100008	M	in: 15" - 17" cm: 38 - 43
100009	L	in: 18" - 20" cm: 46 - 51
100049	XL	in: 21" - 23" cm: 53 - 58



XPF® Caps

Mesh top and wicking material against skin provide maximum comfort for the wearer. Lead equivalency 0.5 mm Pb in front and sides. Machine washable. Ties securely in back. Patented XPF® Technology.

PN	Size	Neck
100098	S	in: 21" - 23" cm: 53 - 58
100099	M	in: 22" - 24" cm: 56 - 61
100100	L	in: 23" - 25" cm: 58 - 64



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com

DISPOSABLE

- 0.5 mm lead equivalency
- Not cleanable, limited use
- Packaged 12/bag



Thyroid Collar

PN	Size	Neck
100001	S	in: 10" - 14" cm: 25 - 35
100002	M	in: 15" - 17" cm: 38 - 43
100003	L	in: 18" - 20" cm: 46 - 51
100051	XL	in: 21" - 23" cm: 53 - 58



Cap

PN	Size	Head
100004	S	in: 21" - 23" cm: 53 - 58
100005	M	in: 22" - 24" cm: 56 - 61
100006	L	in: 23" - 25" cm: 58 - 64