



BLOXR

XPF® Patient Protection

PATIENT SHIELD

Lead equivalency 0.5 mm Pb in front

Designed for placement over radio-sensitive organs and gonads to protect the patient from harmful ionizing scatter radiation.

Patient Shields

PN*	Size	PN*	Size
100201	18" x 18" (46 x 46 cm)	100210	24" x 24" (61 x 61 cm)
100202	20" x 20" (51 x 51 cm)	100211	24" x 36" (61 x 91 cm)
100203	14" x 16" (36 x 41 cm)	100212	24" x 46" (61 x 117 cm)
100204	11" x 34" (28 x 86 cm)	100213	8" x 10" (20 x 25 cm)
100206	16" x 24" (41 x 61 cm)	100214	12" x 12" (31 x 31 cm)
100207	18" x 24" (46 x 61 cm)	100215	16" x 12" (41 x 31 cm)
100208	20" x 24" (51 x 61 cm)	100216	12" x 14" (31 x 36 cm)

Options/Accessories:

Embroidery • Colored binding • Belt with hook/loop closure
Custom sizes available



Orientation is important: the BLOXR® logo on the Patient Shield should face the radiation source

PATIENT PONCHO

Lead equivalency 0.5 mm Pb in front

Minimizes a patient's exposure to ionizing radiation while providing anatomical access for the clinician. The sleeveless design enables it to be donned without moving the patient's arm. Poncho is lightweight so it doesn't put uncomfortable pressure on the patient. Built in flap covers the thyroid area. Extenders can be added for extra length for taller patients. Helpful during orthopedic, interventional, podiatry, dentistry and trauma cases.

Patient Poncho

Extender Sets (optional)

PN*	Size	Size	PN*	Lengths
100231	#1 Small Child	For #1	100240	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm)
100232	#2 Child	For #2	100241	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)
100233	#3 Small Adult	For #3	100242	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)
100234	#4 Adult	For #4	100243	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)
100235	#5 Large Adult	For #5	100244	4" (10.2 cm), 6" (15.2 cm), 8" (20.3 cm), 10" (25.4 cm)

Options/Accessories: Embroidery • Colored binding



FRONT shown with extender



BACK



*Add color code suffix to PN



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com