

Bra Inserts for Breast Protection



BLOXR

Bra Inserts

XPF® Bra Inserts provide 0.5 mm lead equivalent protection over the breast and easily slip inside bra. Useful for those concerned about or at risk for breast cancer. Helpful in situations where the size of the radiation protection apron provided to the user does not adequately cover the chest area.

Lined with soft wicking material for the wearers' comfort. Can be wiped clean or machine washed and dried. Available in standard cup sizes A through D.

XPF Inserts are sold as a pair in hot pink or black.



PN	Size, color
10013A-BLK	Cup size A, black
10013A-PNK	Cup size A, hot pink
10013B-BLK	Cup size B, black
10013B-PNK	Cup size B, hot pink
10013C-BLK	Cup size C, black
10013C-PNK	Cup size C, hot pink
10013D-BLK	Cup size D, black
10013D-PNK	Cup size D, hot pink



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com