

Radiation Protection Eyeglasses



NIKE Brazen

Rx: 0 to +/-4.50 **SV** **BF** **PRO**

PN	Color
7NKBRZ-BLK	Black
7NKBRZ-TOR	Tortoise (shown)

Has rubberized nose bridge

NIKE Flow

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7NKFLW-MANMAR	Matte Night Maroon (shown)
7NKFLW-MAANT	Matte Anthracite
7NKFLW-MATTOR	Matte Tokyo Tortoise

NIKE Legend

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7NKLGD-BLKV	Black Volt
7NKLGD-WGRY	Wolf Gray (shown)
7NKLGD-WHT	White

NIKE Mercurial 6.0

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7NKMEC6.0-NAV	Navy (shown)
7NKMEC6.0-BLK	Black
7NKMEC6.0-GRY	Gray

NIKE Premier 6.0

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7NKPRE6.0-BLK	Black

NIKE Unrest

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7NKUNR-MAANTC	Matte Anthracite Camo (shown)
7NKUNR-MASQBLU	Matte Squadron Blue
7NKUNR-DPWGRY	Deep Pewter Gray

Bolle Keelback

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7BOKE-SBLK	Shiny Black (shown)
7BOKE-SBLKBLU	Shiny Black/Blue Translucent

Bolle King Snake

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7BOKS-TOR	Tortoise (shown)
7BOKS-SBLK	Shiny Black
7BOKS-MABLU	Matte Blue

Bolle Vibe

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7BOVIBE-SBLK	Shiny Black (shown)
7BOVIBE-MAGRY	Matte Gray
7BOVIBE-MABLU	Matte Blue

Smith Optics Frontman

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7SOFM-MABLK	Matte Black

Smith Optics Director

Rx: 0 to +/-3.50 **SV**

PN	Color
7SODIR-BLK	Black

Wiley X P-17

Rx: 0 to +/-2.75 **SV** **BF**

PN	Color
7WYP17-TOR	Tortoise (shown)
7WYP17-BLK	Black

Lenses provide 0.75mm Pb protection. Side shields, if included, provide 0.5mm Pb protection. Lifetime warranty on defects in material and workmanship.

S For smaller size faces **Prescription availability:** **SV** Single vision Rx Add -SV suffix **BF** Bifocal Rx Add -BF suffix **PRO** Progressive Rx Add -PRO suffix



Available options for Nike, Bolle, Smith, Wiley and Oracle frames:
 - Anti-Glare lenses
 - Anti-fog Spray bottle

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
 855-256-9729 . 801-590-9880
 www.bloxr.com . Email: order@bloxr.com

Oracle

Rx: 0 to +/-3.75 **SV** **BF**

PN	Color
7ORA-MABLKS	Matte Black/Silver (shown)
7ORA-MGMR	Metallic Gunmetal/Red

Adjustable rubberized nose bridge

Radiation Protection Eyeglasses



Microlite M8

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7MLM8-MDBLU	Midnight Blue

Microlite MX **S**

Rx: 0 to +/-8.0 **SV** **BF** **PRO**

PN	Color
7MLMX-BLK	Black
7MLMX-NAV	Navy Blue (shown)
7MLMX-CRI	Crimson
7MLMX-SIL	Silver

Microlite MX XL

Rx: 0 to +/-8.0 **SV** **BF** **PRO**

PN	Color
7MLMXL-BLK	Black
7MLMXL-NAV	Navy Blue
7MLMXL-CRI	Crimson (shown)
7MLMXL-SIL	Silver

Microlite Plus **S**

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7MLPLS-BLK	Black
7MLPLS-CRI	Crimson
7MLPLS-PNK	Pink (shown)
7MLPLS-SIL	Silver

Microlite M8 with anti-fog

No Rx

PN	Color
7MLM8-MDBLU-F	Midnight Blue

Comes with Side Shields

Microlite MX with anti-fog **S**

No Rx

PN	Color
7MLMX-BLK-F	Black
7MLMX-NAV-F	Navy Blue (shown)
7MLMX-CRI-F	Crimson
7MLMX-SIL-F	Silver

Comes with Side Shields

Microlite MX XL with anti-fog

No Rx

PN	Color
7MLMXL-BLK-F	Black
7MLMXL-NAV-F	Navy Blue
7MLMXL-CRI-F	Crimson (shown)
7MLMXL-SIL-F	Silver

Comes with Side Shields

Microlite Plus with anti-fog **S**

No Rx

PN	Color
7MLPLS-BLK-F	Black
7MLPLS-CRI-F	Crimson
7MLPLS-PNK-F	Pink (shown)
7MLPLS-SIL-F	Silver

Comes with Side Shields

Lenses provide 0.75mm Pb protection. Side shields, if included, provide 0.5mm Pb protection. Lifetime warranty on defects in material and workmanship.

S For smaller size faces **Prescription availability:** **SV** Single vision Rx Add -SV suffix **BF** Bifocal Rx Add -BF suffix **PRO** Progressive Rx Add -PRO suffix

Fitover style to wear over existing glasses.

Microlite FT

PN	Color
7MLFT-BLK	Black
7MLFT-TOR	Tortoise (shown)

Comes with Side Shields

Microlite FT with anti-fog

PN	Color
7MLFT-BLK-F	Black
7MLFT-TOR-F	Tortoise (shown)

Comes with Side Shields

Large Fitover

PN	Color
7FTO-MABLK	Shiny Black

Comes with Side Shields

Microlite Plus XL

Rx: 0 to +/-3.50 **SV** **BF**

PN	Color
7MLPXL-BLK	Black
7MLPXL-CRI	Crimson
7MLPXL-NAV	Navy Blue (shown)
7MLPXL-SIL	Silver

Microlite Plus XL with anti-fog

No Rx

PN	Color
7MLPXL-BLK-F	Black
7MLPXL-CRI-F	Crimson
7MLPXL-NAV-F	Navy Blue (shown)
7MLPXL-SIL-F	Silver

Comes with Side Shields



BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
 855-256-9729 . 801-590-9880
 www.bloxr.com . Email: order@bloxr.com

Options available for Microlite frames:
 - Name engraving on temple bar
 - Super Hydrophic & Anti-Reflective coating
 - Fog Buster Wipes, box of 60 packets
 Laser eyewear also available.

Microlite is a registered trademark of K&R Medical, LLC