

BLOXR® Reusable Facemask

- Choose from 4 sizes: Adult, Small Adult/Teen, Child (age 7-12), Small Child (age 3-6)
- Outer fabric is nylon ripstop finished with a polyurethane coating and a water repellent finish. The inner fabric is a breathable, moisture-wicking polyester mesh.
- Facemask material tested to be comparable or better than N95 in aerosol transmission.
- No X-Ray attenuation properties.
- Disinfect with disinfectant wipe or soap and water.
- Machine-washable.
- Designed for extended use.

Disclaimer: Non-sterile. This is not a medical device cleared by the FDA. This option is for when Personal Protective Equipment (N95 recommended for Covid-19) is unavailable. This facemask is manufactured per CDC guidelines for homemade masks. Homemade masks are not considered PPE. Wearing this mask over a N95 mask may extend the life of your N95 mask.


You may clean this mask with a disinfectant wipe, by hand with soap and water, or in the washing machine. It is recommended to hang dry (do not expose to high heat).

Ordering Information (qty. 1 each):

- 100150-BLK Adult, black
- 100151-BLK Small Adult/Teen, black
- 100152-BLK Child (age 7-12), black
- 100153-BLK Small Child (age 3-6), black

*Custom logos
can be added.
Please inquire.*

**NOW
Available in
4 sizes**



Ties securely in back

Designed without uncomfortable elastic behind the ears to avoid instances of skin irritation or cutting when worn for long periods of time.

Available on: bloxr.com/shop
For larger orders (12+), call 1-855-256-9729.

