



BLOXR

XPF® Scatter Radiation Protection Pregnancy Options

► Shield with a Strap Belt

Add an additional layer of protection by wearing an XPF Shield either on top of or underneath an existing garment. Our Shield provides 0.5 mm Pb equivalency and comes in a variety of lengths and widths. Belt strap is adjustable.



Optional embroidery available



► Shield attached to an XPF Skirt or Apron

When purchasing an XPF Skirt or Apron, add a detachable Shield to the inside of the garment. Provides an additional 0.5 mm Pb equivalency near the uterus and is easily removable when not needed.



Shield attaches inside and can be removed

Contact your local representative or BLOXR Customer Service for details.



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com