




# Radiation Protection Apparel MEASURING & SIZING CHART

When determining your body's measurements, refer to the diagram below. Select the appropriate style and size from the list of XPF® attenuating garments.



Circumference:

Head: above ears

Neck: around center of neck

Armhole: top of shoulder aligned with armpit, down the front to 1" below armpit, up the back to top of shoulder

Chest/Bust: at largest point

Waist: at belt area

Low Waist: between the belt area and hips, about 4" below waist

Hips: largest point across backside

Length:

Vest: top of shoulder to 4" below waist

Apron: top of shoulder to 1" below center of knee

Skirt: waist to 1" below center of knee

Scan for  
Frontal Apron  
measuring video



Scan for  
Vest/Skirt  
measuring video



## SKIRT

Size	Waist	Hips	Length
XS	26" - 29"	36"	19"
S	29" - 32"	38"	20"
M	32" - 34"	41"	21"
L	34" - 36"	44"	22"
XL	37" - 40"	49"	23"
XXL	40" - 43"	52"	24"

## VEST

Size	Chest	Low Waist	Armhole	Length
XS	30" - 32"	32" - 34"	18"	21"
S	32" - 35"	37" - 39"	20.5"	22"
M	36" - 41"	39" - 42"	22"	24"
L	41" - 46"	42" - 44"	23"	25"
XL	46" - 51"	44" - 47"	24"	27"
XXL	51" - 53"	48" - 51"	26.5"	28"

## FRONTAL APRON

Size	Chest	Armhole	Length	Max Hgt.
XS	30" - 34"	17.25"	32"	5'0"
S	34" - 38"	20.5"	35"	5'5"
M	38" - 42"	25"	40"	5'10"
L	42" - 46"	28"	44"	6'3"
XL	46" - 50"	30"	46"	6'8"
XXL	50" - 56"	32.5"	47"	6'10"

## THYROID COLLAR

Size	Neck
S	10" - 14"
M	15" - 17"
L	18" - 20"
XL	21" - 23"

## CAP

Size	Head
S	21" - 23"
M	22" - 24"
L	23" - 25"

## VEST with ELASTIC BACK

Size	Chest	Low Waist	Armhole	Length
XS	30" - 32"	32" - 34"	17.5"	21"
S	32" - 35"	37" - 39"	19"	22"
M	36" - 41"	39" - 42"	20.5"	24"
L	41" - 46"	42" - 44"	22.5"	25"
XL	46" - 51"	44" - 47"	24"	27"
XXL	51" - 53"	48" - 51"	25.5"	28"

## FRONTAL APRON with ELASTIC BACK

Size	Chest	Armhole	Length	Max Hgt.
XS	30" - 34"	17.5"	32"	5'0"
S	34" - 38"	21.5"	35"	5'5"
M	38" - 42"	25"	40"	5'10"
L	42" - 46"	28"	44"	6'3"
XL	46" - 50"	30.5"	46"	6'8"
XXL	50" - 56"	34"	47"	6'10"

## DROP APRON

Size	Chest	Length	Max Hgt.
S	34" - 38"	37"	5'5"
M	38" - 42"	42"	5'10"
L	42" - 46"	47"	6'3"
XL	46" - 51"	52"	6'8"

## WRAP APRON

Size	Chest	Armhole	Length	Max Hgt.
XS	30" - 32"	20"	32"	5'0"
S	32" - 35"	21"	35"	5'5"
M	36" - 41"	22"	40"	5'10"
L	41" - 46"	23"	44"	6'3"
XL	46" - 51"	24"	46"	6'8"

## FRONTAL APRON with CLOSED ELASTIC BACK

Size	Chest	Armhole	Length	Max Hgt.
XS	30" - 34"	18"	32"	5'0"
S	34" - 38"	19.5"	35"	5'5"
M	38" - 42"	21.5"	40"	5'10"
L	42" - 46"	23"	44"	6'3"
XL	46" - 50"	24.5"	46"	6'8"
XXL	50" - 56"	26"	47"	6'10"

Not sure? We can help you choose the appropriate size or assist with taking measurements. Custom sizing is available by contacting your local representative.

**BLOXR Solutions, LLC**

2250 S. 1300 West, Suite C, Salt Lake City, UT 84119 . 855-256-9729 . Fax 877-254-4888

Email [order@bloxr.com](mailto:order@bloxr.com) . [www.BLOXR.com](http://www.BLOXR.com)

Made in USA  
U.S. Patent No. 8,993,989  
and 8,754,389