



# BLOXR® XPF® Garment Attributes

- ✓ Standard  
✗ Not available/not applicable  
✚ Available Option (can be added to garment)  
Additional charges may apply.

	Vest	Vest with Elastic Back	Skirt	Frontal Apron	Frontal Apron with Elastic Back	Frontal Apron with Closed Elastic Back	Wrap Apron	Drop Apron
Lead Equivalency	✓ 0.5 mm in front; 0.25 mm in back	✓ 0.5 mm in front; 0.25 mm in back	✓ 0.5 mm in front; 0.25 mm in back	✓ 0.5 mm in front	✓ 0.5 mm in front	✓ 0.5 mm in front; 0.25 mm in back	✓ 0.5 mm in front; 0.25 mm in back	✓ 0.5 mm in front
Standard Sizes	✓ XS, S, M, L, XL, XXL	✓ XS, S, M, L, XL, XXL	✓ XS, S, M, L, XL, XXL	✓ XS, S, M, L, XL, XXL	✓ XS, S, M, L, XL, XXL	✓ XS, S, M, L, XL, XXL	✓ XS, S, M, L, XL	✓ S, M, L, XL
Custom Sizing	✚	✚	✚	✚	✚	✚	✚	✚
Adjust Length	✚	✚	✚	✚	✚	✚	✚	✚
Adjust Overlap	✚	✗	✚	✗	✗	✗	✚	✗
Shoulder Pads	✓	✓	✗	✓	✓	✓	✓	✓
Pocket	✓ on left chest	✓ on left chest	✓ on left hip	✓ on left chest and both hips	✓ on left chest	✓ on left chest	✓ on left chest and left hip	✓ on left chest
Tethering (600001)	✚	✚	✚	✚	✚	✚	✚	✚
Direct Embroidery (600004 or 600005)	✚	✚	✚	✚	✚	✚	✚	✚
Patch Embroidery (600002 or 600003)	✚	✚	✚	✚	✚	✚	✚	✚
Embroidered Logo (600008 and 600009)	✚	✚	✚	✚	✚	✚	✚	✚
6" Elastic Support Belt (100304 - 100309)	✚	✗ already has built-in elastic	✗	✚ replaces standard belt	✗ already has built-in elastic	✗ already has built-in elastic	✚ replaces standard belt	✗
Sleeves (SLV-LEFT, SLV-RIGHT, SLV-PAIR, 600015 or 600016)	✚	✚	✗	✚	✚	✚	✚	✗
Alternate Binding Edge Color (600012)	✚	✚	✚	✚	✚	✚	✚	✚
Dosimeter Pocket (600007)	✚	✚	✚	✚	✚	✚	✚	✚
Quick Release Drop Shoulder (600010)	✗	✗	✗	✗	✚ right, left or both shoulders	✗	✗	✗ has a "drop" design
Back Shoulder Strap (600011)	✗	✗	✗	✓ already includes back strap	✚	✗	✗	✗
Anti-Slip Lining (600006)	✗	✗	✚	✗	✗	✗	✗	✗
Inner Lumbar Band (600014)	✗	✗	✚	✗	✗	✗	✚	✗
Suspenders (100114)	✗	✗	✚	✗	✗	✗	✗	✗

	Thyroid Collar	High Neck Thyroid Collar	Max Thyroid Collar	Cap
Lead Equivalency	✓ 0.5 mm in front over thyroid	✓ 0.5 mm in front over full bib and sides	✓ 0.5 mm in front over full bib and sides	✓ 0.5 mm in front and sides
Standard Sizes	✓ S, M, L, XL	✓ S, M, L	✓ one size fits most	✓ S, M, L
Tethering (600001)	✚	✚	✚	✚
Non-standard closure (magnet 600017 or clip 600018)	✚	✗ No magnet (clip only)	✚	✗
Direct Embroidery (600004 or 600005)	✚	✚	✚	✚
Patch Embroidery (600002 or 600003)	✚	✚	✚	✚
Embroidered Logo (600008 and 600009)	✚	✚	✚	✚
Alternate Binding Edge Color (600013)	✗	✚	✗	✚
Dosimeter Pocket (600007)	✚	✚	✚	✚



Scan code for BLOXR sales literature and additional information.



BLOXR Solutions, LLC  
2250 S. 1300 West, Suite C  
Salt Lake City, UT 84119 USA  
855-256-9729 . 801-590-9880  
www.bloxr.com  
Email: order@bloxr.com