

Additional XPF® Items



BLOXR

PROTECTIVE SHIELDS

Lead equivalency 0.5 mm Pb in front

Multi-sized and multi-use shields designed for placement over radio-sensitive organs and gonads for protection from harmful ionizing scatter radiation.

Protective Panel Shields

PN*	Size	PN*	Size
100201	18" x 18" (46 x 46 cm)	100210	24" x 24" (61 x 61 cm)
100202	20" x 20" (51 x 51 cm)	100211	24" x 36" (61 x 91 cm)
100203	14" x 16" (36 x 41 cm)	100212	24" x 46" (61 x 117 cm)
100204	11" x 34" (28 x 86 cm)	100213	8" x 10" (20 x 25 cm)
100206	16" x 24" (41 x 61 cm)	100214	12" x 12" (31 x 31 cm)
100207	18" x 24" (46 x 61 cm)	100215	16" x 12" (41 x 31 cm)
100208	20" x 24" (51 x 61 cm)	100216	12" x 14" (31 x 36 cm)

Options: ✓ Monogram/Embroidery
✓ Colored Binding Edge ✓ Strap belt

Ask about options for Pregnancy Shields.

Orientation is important: the BLOXR® logo on the Protective Shield should face the radiation source.



PATIENT PONCHO

Lead equivalency 0.5 mm Pb in front

Minimizes a patient's exposure while providing anatomical access for the clinician. Sleeveless design enables donning without moving the patient's arm. Has built in thyroid coverage and optional extenders. Helpful for orthopedic, interventional, podiatry, dentistry and trauma cases.



FRONT
shown with extender



BACK

Patient Poncho

Patient Poncho		Extender Sets (optional)		
PN*	Size	Size	PN*	Lengths
100231	#1 Small Child	For #1	100240	4", 6", 8"
100232	#2 Child	For #2	100241	4", 6", 8", 10"
100233	#3 Small Adult	For #3	100242	4", 6", 8", 10"
100234	#4 Adult	For #4	100243	4", 6", 8", 10"
100235	#5 Large Adult	For #5	100244	4", 6", 8", 10"

Options:
✓ Embroidery
✓ Colored Binding Edge



*Add the color code suffix to PN.



BLOXR