



BLOXR

XPF® Scatter Radiation Protection Pregnancy Options

► XPF Skirt or XPF Apron with built-in (removable) Pregnancy Shield

These BLOXR XPF Skirts and Aprons are equipped with a detachable Pregnancy Shield that is attached on the inside, at the front of the garment. In addition to the garment's protection, the Shield provides an additional 0.5 mm Pb equivalency near the uterus. The shield is easily removable when not needed.

XPF Skirt with detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area/0.25 mm Pb back.
If shield is removed, 0.5 mm Pb front/0.25 mm Pb back.

Part Number*	Garment Size (shield W x H)
10P106-	Extra Small (15" x 16")
10P027-	Small (16" x 17")
10P028-	Medium (17" x 18")
10P029-	Large (18" x 20")
10P039-	X-Large (19" x 22")
10P102-	XX-Large (19" x 22")

XPF Frontal Apron with detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area.
If shield is removed, 0.5 mm Pb in front.

Part Number*	Garment Size (shield W x H)
10P019-	Extra Small (17" x 13")
10P010-	Small (19" x 16")
10P011-	Medium (21" x 17")
10P012-	Large (24" x 22")
10P038-	X-Large (26" x 22")
10P045-	XX-Large (28" x 22")

XPF Frontal Apron with Elastic Back and detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area.
If shield is removed, 0.5 mm Pb in front.

Part Number*	Garment Size (shield W x H)
10P420-	Extra Small (17" x 13")
10P421-	Small (19" x 16")
10P422-	Medium (21" x 17")
10P423-	Large (24" x 22")
10P424-	X-Large (26" x 22")
10P425-	XX-Large (28" x 22")

*add color suffix



Contact your local representative or
BLOXR Customer Service for details.



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com