



BLOXR

XPF® Scatter Radiation Protection Pregnancy Options

► Pregnancy Shield with a Strap Belt

Add an additional layer of protection by wearing an XPF Pregnancy Shield over the top of, or underneath, an existing garment. Our Shield provides 0.5 mm Pb equivalency in front. The attached belt strap is made of a webbing material and has velcro-type closures, making it easy to adjust on the body.



Optional
embroidery
available

Pregnancy Shield with attached Strap Belt

Lead equivalency: 0.5 mm Pb in front.

Part Number*	Size (dimensions W x H)
10P290-	Extra Small (17" x 13")
10P291-	Small (19" x 16")
10P292-	Medium (21" x 17")
10P293-	Large (24" x 22")
10P294-	X-Large (26" x 22")
10P295-	XX-Large (28" x 22")

*add color suffix



Color Suffixes: BLU (blue), BLK (black), GRY (gray), GRN (green), RED (red), PUR (purple), TUR (turquoise), CAM (camouflage), ORA (orange), POP (poppy red), YEL (yellow), PNK (hot pink). Colorblock: SKUC (skulls pirate), TROC (tropical), PENC (penguins), PINC (pineapples), FSHC (fish), BEAC (bears), DONC (donuts), HEAC (hearts), LEOC (leopard), TARC (tartan), AVOC (avocado), COOC (cookies), BONC (bones), HOPC (hops), JUNC (jungle flowers), RAIC (rainbow hearts), SPAC (outer space), TIEC (tie dye), SHAC (sharks), USAC (USA stars), CHL (chill), FLAC (flamingo), GEOC (geometric neon), NEMC (clownfish "Nemo"), PIZC (pizza), SSKC (scary skulls), TULC (tulips).



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com