




How to Wear the XPF® Bun Cap

Specifically designed for tucking away long hair.



- 1, 2 - Gather hair into a bun or pony tail.
- 3 - Tuck hair into the extended pouch portion at the rear of the XPF Bun Cap.
- 4 - Grasp the ties to begin securing cap.
- 5 - Tighten the ties.
- 6 - Tie a bow to finish securing the cap on your head.



Scan for video demonstration.



Bun Cap

PN*	Size	Head circumference
100598	S	in: 21"-23" (cm: 53-58)
100599	M	in: 22"-24" (cm: 56-61)
100600	L	in: 23"-25" (cm: 58-64)

*Add color suffix



BLOXR Solutions, LLC
 2250 S. 1300 West, Suite C, Salt Lake City, UT 84119
 801-590-9880 . Toll Free 855-256-9729
 Email: order@bloxr.com | www.bloxr.com

79126A (JAN-20-2026)