

# Radiation Protection Eyeglasses



**BLOXR**

**Prescription availability:**



Single vision Rx  
Add -SV suffix



Bifocal Rx  
Add -BF suffix



Progressive Rx  
Add -PRO suffix



Not Available

**NIKE Brazen**  
Rx: 0 to +/-4.50  SV  BF  PRO

PN: 7NKBRZ-BLK Color: Black (shown)



\*Size: 60-18-134

✓ Rubberized nose bridge



**NIKE Morph**  
Rx: 0 to +/-4.50  SV  BF  PRO

PN: 7NMOR-MBLK Color: Matte Black (shown)



\*Size: 59-16-145

✓ Lower nose bridge (low profile)



**NIKE Premier 6.0**  
Rx: 0 to +/-3.50  SV  BF  PRO

PN: 7NKPRE6.0-BLK Color: Black (shown)




\*Size: 60-13-135



**NIKE 7165**  
Rx: 0 to +/-3.50  SV  BF  PRO

PN: 7NK7165-BLK Color: Black (shown)



\*Size: 52-17-140



**NIKE Radeon 1**  
Rx: 0 to +/-8.00  SV  BF  PRO

PN: 7NKRAD-NAV Color: Navy Blue (shown)  
7NKRAD-MBLK Color: Matte Black



\*Size: 55-17-145



**Bolle King Snake**  
Rx: 0 to +/-3.50  SV  BF  PRO

PN: 7BOKS-TOR Color: Tortoise (shown)  
7BOKS-SBLK Color: Shiny Black



\*Size: 62-17-125




**Wiley X P-17**  
Rx: 0 to +/-2.75  SV  BF  PRO

PN: 7WYP17-TOR Color: Tortoise (shown)  
7WYP17-BLK Color: Black



\*Size: 61-18-120



**Large Fiterover**  
No Rx  SV  BF  PRO

PN: 7FTO-MABLK Color: Shiny Black (shown)



\*Size: 62-15-135  
Frame height: 48 mm

✓ Wear over existing glasses  
✓ Side shields



**Phillips OP23**  
Rx: 0 to +/-7.0  SV  BF  PRO

PN: 7MLOP23-NAV Color: Navy Blue (shown)  
7MLOP23-BLK Color: Black

FF Optional anti-fog lenses (non-Rx only)



\*Size: 54-18-125



**Phillips 15011**  
Rx: 0 to +/-3.00  SV  BF  PRO

PN: 7ML15011-NAV Color: Navy Clear (shown)

FF Optional anti-fog lenses (non-Rx only)



\*Size: 57-15-125

✓ Side shields  
✓ Adjustable nose pads



**Phillips 525**  
Rx: 0 to +/-4.0  SV  BF  PRO

PN: 7ML525-BLK Color: Black (shown)

FF Optional anti-fog lenses (non-Rx only)



\*Size: 59-19-135

✓ Side shields  
✓ Adjustable nose pads  
✓ Metal frame



**Phillips 33 Fiterover**  
No Rx  SV  BF  PRO

PN: 7ML33-TOR Color: Tortoise (shown)  
7ML33-BLK Color: Black

FF Optional anti-fog lenses (non-Rx only)



\*Size: 60-16-12 /  
Frame height: 41 mm

✓ Wear over existing glasses  
✓ Side shields



Anti-Fog (Fog-Free) option is available on Phillips non-prescription (plano) lenses only. Add -FF suffix.

Engraving Option: Personalize with name engraving on the temple bar - add part #7MLENGR (not available on Phillips 525).

**Lenses provide 0.75mm Pb protection. Side shields, if included, provide 0.5mm Pb protection.**

\*Approximate Size (mm) = lens width-bridge-temple length



**BLOXR**

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA  
855-256-9729 . 801-590-9880  
www.bloxr.com . Email: order@bloxr.com

79018L (MAR 30-2026)