



BLOXR

XPF® Scatter Radiation Protection Pregnancy Options

► Pregnancy Shield with a Strap Belt

Add an additional layer of protection by wearing an XPF Pregnancy Shield (Half Apron) over the top of, or underneath, an existing garment. Our Shield provides 0.5 mm Pb equivalency in front. The attached belt strap is made of a webbing material and has velcro-type closures, making it easy to adjust on the body.



Optional
embroidery
available

Pregnancy Shield with attached Strap Belt

Lead equivalency: 0.5 mm Pb in front.

Part Number*	Size (dimensions W x H)
10P290-	Extra Small (17" x 13")
10P291-	Small (19" x 16")
10P292-	Medium (21" x 17")
10P293-	Large (24" x 22")
10P294-	X-Large (26" x 22")
10P295-	XX-Large (28" x 22")

*add color suffix

Color Suffixes: BLU (blue), BLK (black), GRY (gray), GRN (green), RED (red), PUR (purple), TUR (turquoise), CAM (camouflage), ORA (orange), POP (poppy red), YEL (yellow), PNK (hot pink). Colorblock: SKUC (skulls pirate), TROC (tropical), PENC (penguins), PINC (pineapples), FSHC (fish), BEAC (bears), DONC (donuts), HEAC (hearts), LEOC (leopard), TARC (tartan), AVOC (avocado), COOC (cookies), BONC (bones), HOPC (hops), JUNC (jungle flowers), RAIC (rainbow hearts), SPAC (outer space), TIEC (tie dye), SHAC (sharks), USAC (USA stars), CHL (chill), FLAC (flamingo), GEOC (geometric neon), NEMC (clownfish "Nemo"), PIZC (pizza), SSKC (scary skulls), TULC (tulips).



ALSO AVAILABLE IN SETS

Set of 3 Pregnancy Shields with attached Strap Belt
includes embroidery of "Baby on Board" on each shield

Lead equivalency: 0.5 mm Pb in front.

Part Number	Size - Color
10P290SET-AMULTI	1 Small - Hot Pink, 1 Medium - Turquoise, 1 Large - Purple
10P290SET-BMULTI	1 Medium - Yellow, 1 Large - Poppy Red, 1 X-Large - Green



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com



BLOXR

XPF® Scatter Radiation Protection Pregnancy Options

► XPF Skirt or XPF Apron with built-in (removable) Pregnancy Shield

These BLOXR XPF Skirts and Aprons are equipped with a detachable Pregnancy Shield that is attached on the inside, at the front of the garment. In addition to the garment's protection, the Shield provides an additional 0.5 mm Pb equivalency near the uterus. The shield is easily removable when not needed.

XPF Skirt with detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area/0.25 mm Pb back.
If shield is removed, 0.5 mm Pb front/0.25 mm Pb back.

Part Number*	Garment Size (shield W x H)
10P106-	Extra Small (15" x 16")
10P027-	Small (16" x 17")
10P028-	Medium (17" x 18")
10P029-	Large (18" x 20")
10P039-	X-Large (19" x 22")
10P102-	XX-Large (19" x 22")

XPF Frontal Apron with detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area.
If shield is removed, 0.5 mm Pb in front.

Part Number*	Garment Size (shield W x H)
10P019-	Extra Small (17" x 13")
10P010-	Small (19" x 16")
10P011-	Medium (21" x 17")
10P012-	Large (24" x 22")
10P038-	X-Large (26" x 22")
10P045-	XX-Large (28" x 22")

XPF Frontal Apron with Elastic Back and detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area.
If shield is removed, 0.5 mm Pb in front.

Part Number*	Garment Size (shield W x H)
10P420-	Extra Small (17" x 13")
10P421-	Small (19" x 16")
10P422-	Medium (21" x 17")
10P423-	Large (24" x 22")
10P424-	X-Large (26" x 22")
10P425-	XX-Large (28" x 22")

*add color suffix



Contact your local representative or
BLOXR Customer Service for details.



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com