



BLOXR

XPF® Scatter Radiation Protection Maternity Apron

► XPF Maternity Style Frontal Apron with CLOSED Full Length Elastic Back

Provides full front and back coverage.

Has slightly longer length in front for extra protection. Ample velcro ensures adjustability of the straps as the pregnancy term progresses.

Built-in lumbar support and padded shoulders enhance comfort.

Includes one adjustable top shoulder strap and a pocket on the left side of the chest.



► Shown with optional tether. Monogramming available.

Back View



XPF Maternity Style Frontal Apron with CLOSED Full Length Elastic Back

Lead equivalency: 0.5 mm Pb in front and 0.25 mm Pb full back.

Part Number*	Size
10P620-	Extra Small
10P621-	Small
10P622-	Medium
10P623-	Large
10P624-	X-Large
10P625-	XX-Large

*add color suffix

Color Suffixes: BLU (blue), BLK (black), GRY (gray), GRN (green), RED (red), PUR (purple), TUR (turquoise), CAM (camouflage), ORA (orange), POP (poppy red), YEL (yellow), PNK (hot pink). Colorblock: SKUC (skulls pirate), TROC (tropical), PENC (penguins), PINC (pineapples), FSHC (fish), BEAC (bears), DONC (donuts), HEAC (hearts), LEOC (leopard), TARC (tartan), AVOC (avocado), COOC (cookies), BONC (bones), HOPC (hops), JUNC (jungle flowers), RAIC (rainbow hearts), SPAC (outer space), TIEC (tie dye), SHAC (sharks), USAC (USA stars), CHL (chill), FLAC (flamingo), GEOC (geometric neon), NEMC (clownfish "Nemo"), PIZC (pizza), SSKC (scary skulls), TULC (tulips).

Contact your local representative or
BLOXR Customer Service for details.



BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com

BLOXR Solutions is a a certified minority, woman-owned small business.