XPF® Scatter Radiation Protection
Pregnancy Options

Pregnancy Shield with a Strap Belt

Add an additional layer of protection by wearing an XPF Pregnancy Shield (Half Apron) over the top of, or underneath, an existing garment. Our Shield provides 0.5 mm Pb equivalency in front. The attached belt strap is made of a webbing material and has velcro-type closures, making it easy to adjust on the body.





Optional embroidery available





Pregnancy Shield with attached Strap Belt

Lead equivalency: 0.5 mm Pb in front.

Part Number*	Size (dimensions W x H)
10P290-	**Extra Small** (17" x13")
10P291-	**Small** (19" x 16")
10P292-	**Medium** (21" x 17")
10P293-	**Large** (24" x 22")
10P294-	**X-Large** (26" x 22")
10P295-	**XX-Large** (28" x 22")

*add color suffix

Color Suffixes: BLU (blue), BLK (black), GRY (gray), GRN (green), RED (red), PUR (purple), TUR (turquoise), CAM (camouflage, ORA (orange), POP (poppy red), YEL (yellow), PNK (hot pink). Colorblock: SKUC (skulls pirate), TROC (tropical), PENC (penguins), PINC (pineapples), FSHC (fish), BEAC (bears), DONC (donuts), HEAC (hearts), LEOC (leopard), TARC (tartan), AVOC (avocado), COOC (cookies), BONC (bones), HOPC (hops), JUNC (jungle flowers), RAIC (rainbow hearts), SPAC (outer space), TIEC (tie dye), SHAC (sharks), USAC (USA stars), CHL (chill), FLAC (flamingo), GEOC (geometric neon), NEMC (clownfish "Nemo"), PIZC (pizza), SSKC (scary skulls), TULC (tulips).

ALSO AVAILABLE IN SETS

Set of 3 Pregnancy Shields with attached Strap Belt includes embroidery of "Baby on Board" on each shield

Lead equivalency: 0.5 mm Pb in front.

Part Number	Size - Color
10P290SET-AMULTI	**1 Small - Hot Pink, 1 Medium - Turquoise, 1 Large - Purple**
10P290SET-BMULTI	**1 Medium - Yellow, 1 Large - Poppy Red, 1 X-Large - Green**

BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com

XPF® Scatter Radiation Protection
Pregnancy Options

XPF Skirt or XPF Apron with built-in (removable) Pregnancy Shield

These BLOXR XPF Skirts and Aprons are equipped with a detachable Pregnancy Shield that is attached on the inside, at the front of the garment. In addition to the garment's protection, the Shield provides an additional 0.5 mm Pb equivalency near the uterus. The shield is easily removable when not needed.

XPF Skirt with detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area/0.25 mm Pb back.
If shield is removed, 0.5 mm Pb front/0.25 mm Pb back.

Part Number*	Garment Size (shield W x H)
10P106-	**Extra Small** (15" x 16")
10P027-	**Small** (16" x 17")
10P028-	**Medium** (17" x 18")
10P029-	**Large** (18" x 20")
10P039-	**X-Large** (19" x 22")
10P102-	**XX-Large** (19" x 22")

XPF Frontal Apron with detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area.
If shield is removed, 0.5 mm Pb in front.

Part Number*	Garment Size (shield W x H)
10P019-	**Extra Small** (17" x 13")
10P010-	**Small** (19" x 16")
10P011-	**Medium** (21" x 17")
10P012-	**Large** (24" x 22")
10P038-	**X-Large** (26" x 22")
10P045-	**XX-Large** (28" x 22")

XPF Frontal Apron with Elastic Back and detachable Pregnancy Shield

Lead equivalency: 0.5 mm Pb in front with 1.0 mm Pb in shield area.
If shield is removed, 0.5 mm Pb in front.

Part Number*	Garment Size (shield W x H)
10P420-	**Extra Small** (17" x 13")
10P421-	**Small** (19" x 16")
10P422-	**Medium** (21" x 17")
10P423-	**Large** (24" x 22")
10P424-	**X-Large** (26" x 22")
10P425-	**XX-Large** (28" x 22")

*add color suffix





Contact your local representative or BLOXR Customer Service for details.

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com

BLOXR Solutions is a a certified minority, woman-owned small business.

79118D (MAY 22-2026)

XPF® Scatter Radiation Protection
Maternity Apron

XPF Maternity Style Frontal Apron with CLOSED Full Length Elastic Back

Provides full front and back coverage.

Has slightly longer length in front for extra protection. Ample velcro ensures adjustability of the straps as the pregnancy term progresses.

Built-in lumbar support and padded shoulders enhance comfort.

Includes one adjustable top shoulder strap and a pocket on the left side of the chest.



Shown with optional tether. Monogramming available.

Back View



XPF Maternity Style Frontal Apron with CLOSED Full Length Elastic Back

Lead equivalency: 0.5 mm Pb in front and 0.25 mm Pb full back.

Part Number*	Size
10P620-	**Extra Small**
10P621-	**Small**
10P622-	**Medium**
10P623-	**Large**
10P624-	**X-Large**
10P625-	**XX-Large**

*add color suffix

Color Suffixes: BLU (blue), BLK (black), GRY (gray), GRN (green), RED (red), PUR (purple), TUR (turquoise), CAM (camouflage, ORA (orange), POP (poppy red), YEL (yellow), PNK (hot pink). Colorblock: SKUC (skulls pirate), TROC (tropical), PENC (penguins), PINC (pineapples), FSHC (fish), BEAC (bears), DONC (donuts), HEAC (hearts), LEOC (leopard), TARC (tartan), AVOC (avocado), COOC (cookies), BONC (bones), HOPC (hops), JUNC (jungle flowers), RAIC (rainbow hearts), SPAC (outer space), TIEC (tie dye), SHAC (sharks), USAC (USA stars), CHL (chill), FLAC (flamingo), GEOC (geometric neon), NEMC (clownfish "Nemo"), PIZC (pizza), SSKC (scary skulls), TULC (tulips).

Contact your local representative or BLOXR Customer Service for details.

BLOXR

BLOXR Solutions, LLC, P.O. Box 5148, North Branch, NJ 08876 USA
855-256-9729 . 801-590-9880
www.bloxr.com . Email: order@bloxr.com

BLOXR Solutions is a a certified minority, woman-owned small business.

79128A (SEP 01-2026)